



# QFLOW 500™ TITANIUM BOLT- SINGLE LUMEN BOLT KIT

*The optimal solution for single port access*




## Dimensionally Compatible Sensors

QFlow 500™ Perfusion Probe
Depth Electrode(s)
Codman Microsensor ICP
Sophysa ICP (PSO-PTT)

*Disclaimer:* Hemedex makes no claims with respect to the functionality and performance of any sensor other than the Hemedex QFlow 500™ Perfusion Probe.

CE 1639

# QFLOW 500™ TITANIUM BOLT- SINGLE LUMEN BOLT KIT



## Contents / Parts:

1. 2.7mm Drill Bit with Depth Collar
2. Hex Wrench
3. Dura Piercer
4. Cranial Bolt



## FEATURES:

- ➔ Provides single port access
- ➔ Other than a drill, requires no additional tools to insert or remove the bolt
- ➔ Wing nut for easy threading
- ➔ Silicone sealing gasket
- ➔ Made of titanium
- ➔ CT compatible
- ➔ MRI Conditional for static magnetic field of 3.0 T or less

## QFLOW 500™ TITANIUM BOLT KITS

Single Lumen Bolt Kit	Cat. # H0000-3631
Dual Lumen Bolt Kit with Introducer	Cat. # H0000-3642
Quad Lumen Bolt Kit	Cat. # H0000-3644

## QFLOW 500™ PERFUSION PROBE

USA/INTL	Cat. # H0000-1600
----------	-------------------

## BOWMAN PERFUSION MONITOR® \*

MODEL 500	
USA (120V)	Cat. # H0000-0710

Includes : (1 each) Perfusion Monitor, power cord, umbilical cord, printer paper, and instruction manual.

## ACCESSORIES

Probe Dura Piercer & Introducer Kit*	Cat. # H0000-3683
Touhy Borst	Cat. # H0000-3680
Compression Fitting	
Male Luer Cap	Cat. # H0000-3682

\* Indicates Not CE Marked



# HEMEDEX®

564 Main Street, Suite 300  
Waltham, MA 02452

Toll Free: 866.HEMEDEX

P: 617.577.1759

F: 617.577.9328

[salesteam@hemedex.com](mailto:salesteam@hemedex.com)

[www.hemedex.com](http://www.hemedex.com)